

Hardness tester – TH 2050S

24 Positions



Features:

- Advance Micro-controller based user-friendly state-of-the-art product design with alphanumeric splash waterproof polyester soft keys for keyboard. User interactive software in dialogue mode for ease of operation with protection against invalid entries.
- 24 Position Carousel for Automatic measurement of tablet parameters – Width, Thickness, Length/Diameter & Hardness. Fulfills the Pharmacopeial requirement of 20 tablet testing in single run.
- The built in Sweeping blades cleans the tablet path for removal of broken pieces, debris & dust.
- Non-Volatile memory storage of 99 methods with parameters with Alphanumeric entries of Sample Name and Identification Number for authentication.
- Automatic run can be performed continuously up to 1-99 tablets. RUN can be re-started on power failure or on System Reset with last run parameters.
- Linear Displacement 0.8mm/sec with fast advance and return movement of breaking arm saves runtime.
- Daily Auto Incremented Run Number and Factory entered CUSTOMER NAME with Instrument Serial Number on report system makes foolproof and GLP compliant.
- Two Tier – a) ADMIN, and b) USER password protection for method editing and Validation. 25 Users can be set with Name and password.
- Built-in Real Time Clock [RTC] for date and time display and on printout.
- Process Stage and Compression Press no. entry, NMT, NLT and Percentage Limit for detected values of Width, Thickness, Diameter/Length & Hardness with Star mark for out of limit values.
- Validation/Calibration protocol with verification for Hardness – 3 point Load Cell calibration with Standard weights, Width/Length/Diameter and 3 point Thickness Calibration with Standard thickness gauges.
- Quick RUN for Independent Width, Thickness, Diameter/Length, Hardness or All parameter measurement in sequence up to 25 sample tablets with GLP information.
- Error indication on display with Audible alarm helps user to trace the problem.
- Printer Compatibility for Deskjet, Inkjet & Dot-matrix printer.
- Ethernet Connectivity by RE01 Module and EDC80 Data capture Software (Optional).

Reports:

- Selectable Report Format, complying with GLP requirements:
 - a) Report giving Run No., Set parameters and Actual parameters during the test process.
 - b) Report of program parameters for 99 methods, with printout of empty program numbers.
 - c) Validation report for Thickness, Width/Diameter/Length, Tablet Hardness.
- Reports can be obtained even after resetting / power off/ power failure conditions.

LABINDIA®
Where commitment is culture

Specifications:

- **Control:** Micro-controller based [Advanced version of microprocessor]
- **Keyboard:** Alphanumeric splash waterproof polyester soft keys.
- **Display:** 40 x 2 line back lighted liquid crystal display [LCD]
- **Sample carousel – 24 positions**
- **Test Parameter: No. of Tablets = 1 – 99.**
 - 1. Width – Manual Tablet orientation (for tablet length upto 25mm)**
 - Measuring Sensor – Stepper controlled Load Touch sensed
 - Measuring Range – 0.1mm to 20.0mm.
 - Measuring Unit – mm or inch.
 - Measuring Accuracy – ± 0.06 mm.
 - 2. Thickness :**
 - Measuring Sensor – Linear Transducer
 - Measuring Range – 0.1mm to 18.0mm.
 - Measuring Unit – mm or inch.
 - Measuring Accuracy – ± 0.06 mm.
 - 3. Diameter/Length:**
 - Measuring Sensor – Stepper controlled Load Touch sensed
 - Measuring Range – 0.1mm to 32.0mm.
 - Measuring Unit – mm or inch.
 - Measuring Accuracy – ± 0.06 mm.
 - 4. Hardness : Speed Control by Stepper motor**
Linear Displacement 0.8mm/sec with fast advance and return of breaking arm.
 - Measuring Sensor – Load Cell – Strain gauge
 - Measuring Range – 2N to 500 N. (600N – Optional)
 - Measuring Unit – N (Newton's) / Kp (Kilo Ponds) / Kgf Sc (Strong Cobb) / Pound.
 - Measuring Accuracy – ± 1 N.
 - 5. Weight**
 - Manual Entry by User (for reporting purpose)
 - Automatic Balance Interface – Sartorius / Mettler when Automatic Wt. Transfer option is selected.
- **Method Storage:** 99 programs with parameters.
- **Data Storage:** Available with Non-Volatile memory.
- **Report Format:**
 - a) GLP & Pharmacopoeia compliant report with Statistics (Avg, Min, Max and RSD)
 - b) Program parameter Report.
 - c) Validation Reports
- **Output:**
 - a) Centronics parallel port for Dot matrix / Desk jet Printer.
 - b) RS232 C ports for – i) Balance connectivity (Sartorius/Mettler)
ii) PC Connectivity for Data Download.
- **Instrument weight:** 22kgs.
- **Dimensions:** 410W x 280H x 460D mm.
- **Power:** 230Vac, $\pm 10\%$, 50 Hz, 65 Watts. / 110V(60Hz) option available.
- **Environmental Operating Conditions:**
 - a) Operation: Indoor.
 - b) Temperature: 15 to 35°C.
 - c) Humidity: 5 to 90% non-condensing.

LABINDIA®

Labindia Analytical Instruments Pvt. Ltd.

www.labindia.com For Enquiries: sales.mfd@labindia.com

FACTORY: Plot No. EL-72, Electronic Zone, TTC Industrial Area, Thane Belapur Road, Navi Mumbai - 400 705.
Tel.: +91-22-2762 6660/61/62/63/65163480/81 Fax: +91-22-27626664 E-mail: lilplfactory@vsnl.net

APPLICATION & TRAINING CENTRE: R-909, TTC Industrial Area, Thane Belapur Road, Rabale, Navi Mumbai - 400 701. Tel.: +91-22-2760 6955 Fax: +91-22-2760 6706

HEAD OFFICE: Thane: 201, Nand Chambers, LBS Marg, Near Vandana Cinema, Thane (W) - 400 602.
Tel.: +91-22-2598 6075 / 2598 6000 Fax: +91-22-2541 0420 / 2533 5940 E-mail: sales.mfd@labindia.com

REGIONAL OFFICES:

Delhi: +91-11-4330 6001 / 10

Chennai: Tel.: +91-44-2434 7008 / 2432 0352

Kolkatta: Tel.: +91-33-2466 1396 / 2466 3362

BRANCH OFFICES:

Bangalore: Tel.: +91-80-2323 0923 / 21

Lucknow: Tel.: +91-522-2346 535 / 2346 496

Hyderabad: Tel.: +91-40-6550 4133 / 2331 3300

Gurgaon: Tel.: +91-124-2843 300 / 2843 600

Tiruvanthapuram: Tel.: +91-471-2324 064 / 2320 082

Vadodara: Tel.: +91-265-2342 929

Pune: Tel.: 098220 43622

Chandigarh: Tel.: +91-172-4090 001-07 / 4090 009



Striving to become the best individuals, we endeavour to foster the best team. Performing sensibly, we try to achieve the best efficiency. Working innovatively, we seek to make the best products. Listening patiently, we excel to offer the best service. So, no matter what your needs are, come to us, GET THE BEST

LABINDIA reserves the right to change specification without notice as part of its continuous programme of product development.