

VACUUM LEAK TESTER

LT 1001



LABINDIA Vacuum Leak Tester 1001 is designed with utmost precision using top notch quality material and advanced equipments. Our Leak Tester comprises a polycarbonate desiccators housing to sustain the vacuum for long time. Leak Tester apparatus finds wide application in numerous industries such as food, pharmaceuticals, chemicals etc. It's also used to test the quality of packing processes in strips, blisters and sachets containing tablets, granulates and liquids.

Salient Features :

- 128 X 64 Monochrome Graphical Display
- Alphanumeric Keypad
- 20 Test and 10 Methods with date and time available.
- Supplied in SS body
- Selectable Pressure Units
- LAN, RS - 232, Parallel Printer, External 80/48 Column and 40 Column Built in Mini Dot Matrix Printer available

LABINDIA
ANALYTICAL
WHERE COMMITMENT IS CULTURE

SPECIFICATIONS

Vacuum Range	Upto 600 mm Hg
Vacuum Setting	Better Vacuum Stabilization Set Level
Pressure Unit	mm Hg / Psi / kPa
Desiccator	8 inch , 10 inch and 12 inch
Display	128 X 64 Monochrome Graphic Display
Keypad	Alphanumeric Keypad
Interface Port	LAN , RS - 232 ,Parallel Printer , External 80/48 Column and 40 Column Built in Mini Dot Matrix Printer
Hold Time	Upto 30 minutes
Power	AC 230V \pm 10V , 50Hz
Accessories	Polycarbonate desiccator along with connecting tube. Instruction Manual
Date Format	USA / ASIA
Real Time Clock	Yes
Calibration	Easy Calibration
Validation	Easy Validation
Batch Test Results	20 Test and 10 Methods with date and time available
Print Report	Standard / User Defined
Password Protection	3 Tier Password Protection
Self Test	Facility to Scan (Valves , Pump and Sensor)
Balance Interface	At Extra Cost (Optional)

*Specifications are subject to change without prior notice

LABINDIA

ANALYTICAL
WHERE COMMITMENT IS CULTURE

www.labindia-analytical.com / For Enquiries: sales.mfd@labindia.com

Labindia Analytical Instruments Pvt. Ltd.

HEAD OFFICE:

Thane: 201, Nand Chambers, LBS Marg, Near Vandana Cinema, Thane (W) - 400 602.
Tel.: +91-22-2598 6000 / 2598 6110 Fax: +91-22-2541 0420 / 2533 5940 E-mail: sales.mfd@labindia.com

REGIONAL OFFICES:

Delhi: Tel.: +91-11-4330 6001 / 10

Chennai: Tel.: +91-44-2434 7008 / 2432 0352

Kolkatta: Tel.: +91-33-2466 3362

BRANCH OFFICES:

Bangalore: Tel.: +91-80-2323 0919 / 20 / 21

Lucknow: Tel.: +91-522-2346 535 / 2346 496

Hyderabad: Tel.: +91-40-4444 3456

Tiruvanthapuram: Tel.: +91-471-2324 064 / 2320 082

Vadodara: Tel.: +91-265-2986 005 / 2986 006

Pune: Tel.: +91-20-2538 0090

Chandigarh: Tel.: +91-172-4090 001-07 / 4090 009

FACTORY: Plot No. EL-72, Electronic Zone, TTC Industrial Area, Thane Belapur Road, Navi Mumbai - 400 705.

Tel.: +91-22-2762 6660 / 61 / 62 / 63 / 6516 3480 / 81 / 6107 6666 Fax: +91-22-27626664 E-mail: lipfactory@vsnl.net

APPLICATION & TRAINING CENTRE: R-909, TTC Industrial Area, Thane Belapur Road, Rabale, Navi Mumbai - 400 701.

Tel.: +91-22-2760 6955 / 6107 6555 Fax: +91-22-2760 6706



Striving to become the best individuals, we endeavour to foster the best team. Performing sensibly, we try to achieve the best efficiency. Working innovatively, we seek to make the best products. Listening patiently, we excel to offer the best service. So, no matter what you needs are, come to us. **GET THE BEST**

LABINDIA reserve the right to change specification without notice as part of its continuous programme of product development.